CC DISK NF CoCr

SI Navodila za uporabo / GB Directions for use
D Gebrauchsanweisung / CZ Návod k použití
HR Uputa za uporabu

CZ
INTERDENT s.r.o.
Foerstrova 12, Strašnice
CZ-10000 Praha
tel. +420/274 814 404
fax +420/274 820 130

HR
INTERDENT d.o.o.
Vinogradski odvojak 2d
HR-10431 Sveta Nedelja
tel. +385 1 3873 644
fax: +385 1 3873 617
mail: interdent@interdent.hr

SRB
INTERDENT d.o.o.
Nehruova 44
SRB-11070 Novi Beograd
tel. +381/11 217 53 74
fax +381/11 3 01 75 33

MK
INTERDENT d.o.o.
Bojmija 1-12M
MK-1000 Skopje
tel. +389/224 628 55
fax +389/224 634 12

SK
INTERDENT SK, s.r.o.
Liščie údolie 57
SK-84231 Bratislava
tel.č. 00421 2 5440554,
00421 2 5440555
e-mail: interdent@interdent.sk

Proizvajalec/Producer/Hersteller
SI - 3000 CELJE · Opekarniška cesta 26
Tel.: +386 (0)3 425-62-00 · Fax: +386 (0)3 425-62-10
E-mail: info@interdent.cc · www.interdent.cc

CE 0197

06/2015
Made in Slovenia

SI

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GB

CC DISK NF CoCr is a casted plate intended for making prosthetic substitutes in CAD/CAM milling machines. CC DISK NF CoCr is a dental alloy without nickel, beryllium, indium and gallium and fulfills the recommendation of the standard EN ISO 22674 and EN ISO 9693-1 for dental alloys. Vickers hardness 285 (HV10) value leads to good milling and drilling behavior and gives perfect polished surface. Only few oxides are formed during ceramic firing. CC DISK NF CoCr is a medical device according European directive 93/42/EEC and is marked by the producer with CE 0197.

COMPOSITION (mass %)		TECHNICAL DATA (Orientation values)	
Co	63	Type	4
Cr	24	Density	8,3 g/cm ³
W	8	Vickers hardness HV 10	285 (HV 10)
Mo	3	Coefficient of thermal expansion	25 - 500°C 13,9 x 10 ⁻⁶ K ⁻¹ 20 - 600°C 14,0x 10 ⁻⁶ K ⁻¹
Si	1	0,2%-Elongation limit	490 MPa (N/mm ²)
Nb	<1	E-modul	ca. 210000 MPa (N/mm ²)
		Ductile yield A5	10%

Ceramic firing

You can use all types of ceramic which you also use in conventional method on cast crowns and bridges (VITA, Ivoclar, Ducera, etc.).

Please be careful about cooling procedure, Co- based alloys must be cooled down slowly after each firing step.

1. Remove connectors and smooth sharp edges. Sand blast the surface with 150 or 250 microns aluminum oxide (Interlox 250 micr or Interlox 150 micr).

2. Odlitek očistite v destilirani vodi s pomočjo ultrazvoka, ter ga razmastite z etilacetatom ali s paro.
3. Oksidacijska peka (samo za kontrolo površine odlitka) naj poteka **5 min. pri cca. 960°C. Po oksidacijski peki obvezno oksidacijsko površino dobro speskatje** z aluminijevim oksidom 250 ali 150 mikronskim in še enkrat razmastite.
4. Prvi sloj, opaker, nanesite tanko in enakomerno po celi konstrukciji, nato sledite navodilom proizvajalca porcelana.
5. Drugi sloj opakra nanesite tako, da pokrijete kompletno kovino, da ne preseva čez opaker.
6. **Po vsaki stopnji peke ohlajajte počasi.**

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2. Clean the surface in distilled water with ultrasound, degrease with ethyl acetate or under vapor.
3. Carry out oxide firing (only to control the surface) 5 min at 960°C. After oxide firing removing of oxides is required. Use aluminum oxide 250 or 150 micr and again degrease the surface.
4. First layer, opaque, applied thinly and evenly over the construction, then follow direction for use of ceramic producer.
5. Second layer of the opaque apply the way that all the surface of construction is covered and that the metal shade is not visible over the ceramic.
6. **After each firing carry out slowly cooling down.**

Finishing

After last firing step finish visible surface of the alloy with polishing rubbers and polishing paste for Co-Cr alloys (REF 0460 Universal polishing paste for Co-Cr) to high gloss. Remove oxides with glass pearls or very carefully with aluminum oxide 50 microns on the cervical parts of the crowns (that gingival fit is not damaged).

Safety precaution

Wear protective clothing, respiratory mask with filter FFP2, gloves, goggles and use suction units when finishing and sand-blasting the construction.

Warranty!

Our recommendations for use are based upon our own experience and trials and can only be considered a standard values. Our products are subject to further development. Therefore alterations in construction and composition are reserved.

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