

NEW
Milling Machine

All-in-One Titanium Disk Milling

CRAFT PLUS

Servo Motor + Powerful Spindle = Ideal for Metal Disk Milling

Key Features

All-in-One



All-in-one system,
now with titanium
milling performance

Servo Motor & Powerful Spindle



High-precision,
high-speed
metal milling

35° Milling



From full-arch
to complex cases with
exceptional accuracy

Milling Materials



Zirconia

Ceramic

Pre-milled

Titan Disk

Specification

Material	Zirconia, PMMA, G/C, H/C, Premill, Titan Disk
Spindle	SycoTec DC
Spindle Power	Max. 1.37kW / 60,000 RPM
Tool Quantity	21ea
Dimensions (WxDxH)	Machine 545 × 670 × 670mm Table 545 × 670 × 970mm
Weight	Total 145kg
CAM	hyperDENT



DOF Inc. #601-603, 77, Seongsuil-ro, Seongdong-gu, Seoul, 04790 Korea (070-5057-3518)

Copyright © 2025 DOF Inc. All rights reserved.

Design and specifications are subject to change without notice to improve performance.

All-in-One Titanium Disc Milling

CRAFT PLUS