

Ultra-high-speed Scan

Up3d Scanner Up400



Specializing in developing scanners over 8 years

UP400Ultra-high-speed Scan

Automatic add-scan





Open System

Support PLY format Enables to import 3shape and EXOCAD files directly



High Speed Scan

Only need 14s to scan unsectioned model



Upgraded Upscan Software

Optimize the UI layout and upgrade the core algorithm to delive faster scan speed and less post-processing time.

Ultra-high-speed Model Scanning 14s For Unsectioned Model, 11s For Sectioned Model

Powered by upgraded software algorithm, UP400 is able to scan the modle with high speed and improve your efficiency



Unsectioned model




UP400 scan data

UP3D Scanner UP400

Open System, Better Integration Highly Restore The True Color Of The Model To Help Rpd Design

UP400 is able to provide you the ture color images, Combining with variety of design software such as exocad and 3shape to achieve the easiest PRD design





With The Ultra-high-speed Scan Up400 Provides You The Better Scan Experience

Post & core impression scanning

Only need 33 seconds





Triple Tray Impressions
Scanning Only
Need 70 Seconds





UP3D Scanner UPSCAN

UPSCAN

Upgrade Scan Software

New UI interface, post-processing speed upgrade





Freely Switch Scanning Steps

Support freely switch scanning steps, no need to re-create the order, simple and convenient operation



Fill Hole On Abutment

No need to add the wax manually, the software fills the screw holes automatically



Full Arch Automatic Add-scan

One click on the software can add-scan automatically



Mark Margin Line

Mark the margin line directly in the software after scanning, reducing time of the tooth design.

Up3d Scanner Series Parameters

Model UP400

Size 355*268*384 (mm)

Weight 10kg

Cameras 2*1.3 (mp) Accuracy 8 microns

Projector Customized blue-light /

Multi-color scanning

Interface USB 3.0

Output Formats STL; UM; PLY; OBJ Power Supply AC110-240V, 50HZ

Product Features All-in-one scan

One click scanning triple tray

Detect margin line

True color texture scan

Main Feature Unsectioned models scan

All-in-one

smart scan

Scan Speed Full arch: 11s

Unsectioned models scan: 14s

post&core scan: 33s

All-in-one: 20s Triple tray: 70s

True color texture scan: 16s





Digital technology improves life



www.up3ds.com info@up3d.cn +86-755-26983202

ADD: 17915 Sky Park Circle, Bldg.27, Suite E, Irvine, CA 92614, USA