



UP600

Intraoral Scanner

No need spray

Clear margin line, high precision of real color

Intelligent scanning

Smart heating and anti-fogging

Easy learn and use, simple operation

STL, PLY, OBJ open data format

Intelligent hole-filling, with more accurate data

Multiple scanning head to satisfy
different kinds of complicated cases

Impression scanning



Comfortable, accurate, efficient and vivid, building a precision era of digital dentistry

Intraoral Scanner Up600

Combined intelligently with implantodontics, orthodontics and prosthodontics to replace the traditional impression taking

Upgrade to digital dentistry completely with higher accuracy



Shortcut To Digital Treatment

The best option to enter digital dentistry era with the latest IOS scanning technology

A surprise of pleasant starts from a flash of initiation!"

Detailed configuration	UP600
Dimension	275*49.5*47mm
Weight	450g
Depth Of Field	11mm
Scanning Height	0-17mm
Scanning Range	12*11mm
Accuracy	20μm
Scanning Speed	18-20 PFS
Scanning Head Anti-fogging Method	Hot wind heating
Scanning Head	45°, 17.4mm height
	35°, 17.6mm height

Digital technology improves life

ADD: 17915 Sky Park Circle, Bldg.27, Suite E, Irvine, CA 92614, USA

