


# Smart Operating System

## The Dental Sintering Furnace, F20




### The Dental Sintering Furnace, F20

The easy entry to the denture sintering system



#### Power Failure Renewal

Continued sintering after abnormal power failure ensures efficient sintering



#### Multiple Sintering Processes

Professional 20 sintering programs are available



#### Large Capacity

Support 150 crowns sintering simultaneously on three sintering hold racks.



# F20 Sintering Furnace

## Smart Operating System

Model	F20
Dimension	512X352X656(mm)
Weight	64Kg
Rated Voltage	220V±10% 110v±10%
Rated Frequency	50HZ
Fuse	2X~500V 32A
Rated Power	2000W
Heating Element	Four Ultrapure Silicon Molybdenum Bar
Operation Mode	7 - Inch Color Touch Screen
Maximum Temperature	1600°C
Temperature Sensor	Model-B Thermocouple



[www.up3ds.com](http://www.up3ds.com)

[info@up3d.cn](mailto:info@up3d.cn)

+86-755-26983202

ADD: 17915 Sky Park Circle, Bldg.27,  
Suite E, Irvine, CA 92614, USA

